



Mallet Splint from Orficast!

*Instructions: Fold a 6" piece of Orficast in half.
Heat to stick sides together and activate.
Cut lengthwise in middle, leaving lower 1" intact.
Place this uncut end over middle phalanx.
Stretch narrow strips around and under distal phalanx.
Support distal phalanx in slight hyper-extension while
pinching Orficast strips together.
Trim corners and apply Velcro hook with heat gun to
Orficast. Finish with strap over middle phalanx.*



*Deborah A. Schwartz, OTD, OTR/L, CHT
Product and Educational Specialist
Orfit Industries America
Debby.schwartz@orfit.com*

*Advantages:
Quick Mallet Splint!
Easy to do!!
Use Orficast!
Check out www.orfit.com !*