

# Phyaction I

*Phyaction I, for electrotherapy, only with tens -based and interferential currents*

With the Phyaction series you treat your patients in a simple and effective way.

- Large and crystal-clear LCD screen
- Knob for precise intensity adjustment
- Simple and easy to use
- Extensive indication menu
- Compact unit and unique design



Article number: 323 265

## Characteristics

- Therapies :
  - Electro therapy, 2 and 4 poles, 2 channels completely independent
  - Simultaneous therapy (2 different indications treated simultaneously by using electro)
- Currents :
  - Current forms: see on the next page
- Functionalities :
  - Objectives: 0
  - Indication: 84
  - Diagnostic: 0
  - Free memory: 50
- Phyaction I can be connected to Phyaction V

## Technical specifications

Dimensions Phyaction	: 265 x 275 x 122 mm
Weight Phyaction	: 3.650 kg
Weight including accessories	: 4.6 kg
Mains voltage	: 100 - 240 VAC, 50-60 Hz
Maximum power, in operation	: 85 VA
Safety class	: Class I (earthed socket required)
Insulation	: Type BF (floating patient circuit)
Fuses	: 2 x T2AL250 V
Treatment time	: 0 - 60 min.
Languages	: English, Dutch, German, French, Italian, Spanish and Portuguese

## Standard accessories

108 725	Electrode cables
109 959	Electrode rubber 6 x 8 cm
100 658	Chamex bag for electrode 6 x 8 cm (per 4)
108 935	Fixing strap, elastic 5 x 60 cm
100 689	Mains lead
108 919	Test plug F/F 4 mm
115 684	VAS

## Manual

322 857	English
323 000	Safety instruction

## Optional accessories

325 068	BAG (carrying case)
107 348	Vaginal probe
107 349	Anal probe
112 166	Rectal probe
109 958	Electrode rubber 4 x 6 cm (per 2)
109 960	Electrode rubber 8 x 12 cm (per 2)
100 657	Chamex bag for electrode 4 x 6 cm (per 4)
100 659	Chamex bag for electrode 8 x 12 cm (per 4)
326 799	Electrode, adhesive, reusable - ø 3 cm (per 4)
326 810	Electrode, adhesive, reusable - 2,5 x 5 cm (per 4)
326 821	Electrode, adhesive, reusable - 5 x 5 cm (per 4)
326 832	Electrode, adhesive, reusable - 5 x 10 cm (per 4)
328 504	Adhesive electrode kit
113 334	Adapter for adhesive electrode 4 > 2 mm
114 142	Pen electrode with chamex bag
109 944	chamex for pen electrode 10 pieces
108 934	Fixing strap 5 x 30 cm
108 936	Fixing strap 5 x 120 cm
117 210	Mobile 700

# Phyaction I : Current forms



## Tens Currents

<p>Conventional Tens example: Symmetrical</p>	<p>Low Frequency Tens example: Symmetrical</p>	<p>Random frequency Tens example: symmetrical</p>
<p>Burst Tens example: symmetrical</p>		

## NMES Currents

<p>Biphasic Surge example: a-symmetrical</p>	<p>Biphasic Surge example: symmetrical</p>	<p>Intrapulse Interval Surge Symmetrical</p>
<p>2-pole Medium Frequency Surge</p>	<p>Russian Stimulation</p>	<p>Isoplanar Vector field surge</p>

## Interferential Currents

<p>2-pole Medium frequency</p>	<p>Dipole vectorfield</p>	<p>Isoplanar Vectorfield</p>