





SPIO®

Sizing Guidelines

Measurement Key



SPIO® Lower Body Orthosis

Waist circumference is measured at belly button.

Length is from waist to Ankle (P) or bottom of Patella (SP).

Short #/Long #

SP41 / P41	Waist Circumference: 41 cm, SP27.25 cm long / P 38.25 cm long
SP43 / P43	Waist Circumference: 43 cm, SP 29.5 cm long / P 42 cm long
SP45 / P45	Waist Circumference: 45 cm, SP31.75 cm long / P45.75 cm long
SP47 / P47	Waist Circumference: 47 cm, SP 34 cm long/ P 49.5 cm long
SP49 / P49	Waist Circumference: 49 cm, SP 36.25 cm / P 53.25 cm long
SP51 / P51	Waist Circumference: 51 cm, SP 36.25 cm / P 55 cm long
SP53 / P53	Waist Circumference: 53 cm, SP 38.75 cm / P 59.5 cm long
SP55 / P55	Waist Circumference: 55 cm, SP 41.25 cm / P 64 cm long
SP57 / P57	Waist Circumference: 57 cm, SP 43.75 cm / P 68.5 cm long
SP59 / P59	Waist Circumference: 59 cm, SP 46.25 cm long / P 73 cm long
SP61 / P61	Waist Circumference: 61 cm, SP 48.75 cm long / P 77.5 cm long
SP63 / P63	Waist Circumference: 63 cm, SP 51.25 cm long / P 79.5 cm long
SP65 / P65	Waist Circumference: 65 cm, SP 53.25 cm long / P 82.5 cm long
SP67 / P67	Waist Circumference: 67 cm, SP 55.25 cm long / P 85.5 cm long
SP69 / P69	Waist Circumference: 69 cm, SP 57.25 cm long / P 88.5 cm long
SP71 / P71	Waist Circumference: 71 cm, SP 59.25 cm long / P 91.5 cm long
SP73 / P73	Waist Circumference: 73 cm, SP 61.25 cm long / P 94.5 cm long

- Chest measurements are taken with patient lying down
- Pant measurements are taken with patient standing
- If in-between sizes refer to length measurements.

SPIO® TLSO

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Item #

V41	Chest Circumference: 41 cm, 27 cm long
V44	Chest Circumference: 44 cm, 30.75 cm long
V47	Chest Circumference: 47 cm, 34.5 cm long
V51	Chest Circumference: 51 cm, 38.25 cm long
V54	Chest Circumference: 54 cm, 44 cm long
V57	Chest Circumference: 57 cm, 47.75 cm long
V60	Chest Circumference: 60 cm, 51.5 cm long
V64	Chest Circumference: 64 cm, 55.25 cm long
V67	Chest Circumference: 67 cm, 59 cm long
V70	Chest Circumference: 70 cm, 62.75 cm long
V73	Chest Circumference: 73 cm, 66.5 cm long
V77	Chest Circumference: 77 cm, 70.25 cm long

SPIO® LSO

Chest circumference is measured at nipple line.
Length is bottom of Rib Cage to Greater Trochanter.

Item #

AB51	Chest Circumference: 51 cm, 15.5 cm long
AB54	Chest Circumference: 54 cm, 17.75 cm long
AB57	Chest Circumference: 57 cm, 21 cm long
AB60	Chest Circumference: 60 cm, 22.25 cm long
AB64	Chest Circumference: 64 cm, 25.5 cm long
AB67	Chest Circumference: 67 cm, 28 cm long
AB70	Chest Circumference: 70 cm, 32 cm long
AB73	Chest Circumference: 73 cm, 34.5 cm long
AB77	Chest Circumference: 77 cm, 37.5 cm long

SPIO® Upper Body Orthosis

Chest circumference is measured at nipple line.
Length is from Acromion to Greater Trochanter.

Short #/Long #

SS41/ S41	Chest Circumference: 41 cm, 29.25 cm long
SS44/ S44	Chest Circumference: 44 cm, 33.25 cm long
SS47/ S47	Chest Circumference: 47 cm, 37.25 cm long
SS51/ S51	Chest Circumference: 51 cm, 41.25 cm long
SS54/ S54	Chest Circumference: 54 cm, 45.25 cm long
SS57 / S57	Chest Circumference: 57 cm, 49.25 cm long
SS60/ S60	Chest Circumference: 60 cm, 53.25 cm long
SS64/ S64	Chest Circumference: 64 cm, 57.25 cm long
SS67/ S67	Chest Circumference: 67 cm, 61.25 cm long
SS70/ S70	Chest Circumference: 70 cm, 65.25 cm long
SS73/ S73	Chest Circumference: 73 cm, 69.25 cm long
SS77/ S77	Chest Circumference: 77 cm, 73.25 cm long

SPIO® Arm Orthosis

Arm circumference is measured at arm pit.

Item #

AS13	Arm Circumference: 13 cm, 19.25 cm long
AS14	Arm Circumference: 14 cm, 21.25 cm long
AS16	Arm Circumference: 16 cm, 23 cm long
AS17	Arm Circumference: 17 cm, 25.5 cm long
AS18	Arm Circumference: 18 cm, 28.25 cm long
AS20	Arm Circumference: 20 cm, 30.75 cm long
AS21	Arm Circumference: 21 cm, 33.25 cm long
AS22	Arm Circumference: 22 cm, 35.75 cm long
AS24	Arm Circumference: 24 cm, 38.5 cm long
AS25	Arm Circumference: 25 cm, 41 cm long
AS26	Arm Circumference: 26 cm, 43.5 cm long
AS28	Arm Circumference: 28 cm, 46 cm long

SPIO® Wrist-Hand Orthosis

Palm width is measured at widest part above thumb.

Item #

G4	Width of Palm: 4 cm, 9.25 cm long
G4.5	Width of Palm: 4.5 cm, 10.5 cm long
G5	Width of Palm: 5 cm, 11 cm long
G5.5	Width of Palm: 5.5 cm, 12 cm long
G6	Width of Palm: 6 cm, 13.25 cm long
G6.5	Width of Palm: 6.5 cm, 14 cm long
G7	Width of Palm: 7 cm, 15.5 cm long
G7.5	Width of Palm: 7.5 cm, 16.25 cm long
G8	Width of Palm: 8 cm, 17.5 cm long

SPIO® Leg Orthosis

Leg circumference is measured at top of thigh.

Item #

LS14	Leg Circumference: 14 cm, 21.75 cm long
LS15	Leg Circumference: 15 cm, 24.5 cm long
LS16	Leg Circumference: 16 cm, 27.25 cm long
LS17	Leg Circumference: 17 cm, 30 cm long
LS18	Leg Circumference: 18 cm, 32.75 cm long
LS19	Leg Circumference: 19 cm, 40 cm long
LS20	Leg Circumference: 20 cm, 43.5 cm long
LS21	Leg Circumference: 21 cm, 46.75 cm long
LS22	Leg Circumference: 22 cm, 54.5 cm long
LS23	Leg Circumference: 23 cm, 56.5 cm long
LS24	Leg Circumference: 24 cm, 58.25 cm long
LS25	Leg Circumference: 25 cm, 60.25 cm long
LS26	Leg Circumference: 26 cm, 62 cm long
LS27	Leg Circumference: 27 cm, 63.75 cm long
LS28	Leg Circumference: 28 cm, 65.75 cm long